

VINYL EXAMINATION GLOVES POWDER FREE

DISRIPTION

Vinyl Gloves are a reliable and versatile choice for various applications that demand hand protection. These gloves are made from natural clear polyvinylchloride (PVC) material, ensuring durability and flexibility. They are an excellent option for individuals with latex allergies, as they are completely latex-free.

1. **Material:** Our Vinyl Gloves are crafted from high-quality polyvinylchloride (PVC) material, which is known for its strength and flexibility. This ensures that the gloves offer reliable protection for your hands.
2. **Colour:** Vinyl examination gloves commonly come in a clear or transparent color.
3. **Powder Free** Powder-free gloves minimize the risk of powder particles contaminating food or surfaces. This is crucial in maintaining hygiene standards in the food industry.
4. **FDA Compliance:** Our Vinyl Gloves adhere to the strict regulations set by the FDA. They fully comply with FDA rule 177.1950 for food contact, making them suitable for use in food-related applications. Additionally, they meet FDA rule 110.10 GMP for food handling, ensuring the highest standards of safety and hygiene.
5. **DINP Manufacturing:** These gloves are manufactured using Di-isononyl phthalate (DINP), which is a widely accepted plasticizer. DINP contributes to the flexibility and durability of the gloves, making them ideal for a range of tasks.
6. **DOP/DEHP-Free:** The Vinyl Gloves are entirely free of DOP (Di-2-ethylhexyl phthalate) and DEHP (Di(2-ethylhexyl) phthalate), which are substances that may be of concern. This ensures that the gloves are safe for use and in compliance with relevant regulations.

Vinyl Gloves are a dependable choice for a variety of industries and applications, including food service, healthcare, laboratory work, cleaning, and more. Their latex-free nature and compliance with FDA guidelines make them a practical and safe option for those in need of protective handwear. Whether you are a healthcare professional, a food handler, or anyone requiring hand protection, our Vinyl Gloves offer the quality and performance you can rely on.

DESIGN & FEATURES:

- **Ambidextrous Design:** Both gloves in a pair are identical, allowing for convenient use on either hand. This design feature ensures versatility and ease of use.
- **Rolled Cuff:** Many vinyl gloves have a rolled or beaded cuff, providing added strength to the gloves and making them easier to put on and take off without rolling back.
- **Smooth Surface:** Vinyl gloves typically have a smooth surface that allows for easy donning and removal and ensures a comfortable fit during extended wear.

- Textured Fingertips (Optional): Some vinyl gloves come with textured fingertips to enhance grip and handling of objects, especially useful in tasks that require precision or grip in slippery conditions.
- Powder-Free Design: Most vinyl gloves are powder-free to reduce the risk of allergic reactions and minimize contamination in sensitive environments. The absence of powder also helps in maintaining a clean working environment.

Sizes and Thickness:

- Sizes: Available in various sizes such as small, medium, large, and extra-large to accommodate different hand sizes for a proper and comfortable fit.
- Thickness: Vinyl gloves come in different thicknesses, typically measured in mils (1 mil = 0.001 inches). Thicker gloves offer more durability and resistance, while thinner gloves provide better tactile sensitivity.

Latex-Free and Allergen-Free:

- Latex-Free: Vinyl gloves are an excellent alternative for individuals with latex allergies or sensitivities, offering a hypoallergenic option for hand protection.

Applications:

- Wide Range of Applications: Vinyl gloves find utility in various industries such as healthcare, food service, laboratory work, janitorial services, beauty and salon industries, cleaning tasks, and more due to their versatility and cost-effectiveness.

Compliance and Certification:

- Regulatory Compliance: These gloves adhere to industry and regulatory standards, including FDA compliance for food contact and relevant safety standards, ensuring their suitability for specific applications.

These design and feature elements collectively contribute to vinyl gloves' effectiveness in providing protection, comfort, and versatility across different professional settings where hand protection is essential.